

영문규격서

COMMODITY DESCRIPTION

관세분류번호 HSK No.	품목번호 Item No.	품명 및 규격 Description	단위 Unit
A79160.06		Ultrasonic Aspirator system	1 Set

A. FEATURE

1. The first Ultrasonic aspirator to effectively address soft tissue and fine bone dissection by coupling longitudinal oscillation and torsional vibration
2. LT technology makes it possible to emulsify bone with no grabbing or kicking due to the non-rotational reciprocating motion, and oscillation.
3. Frangible tip covers for improved visualization and console automatically recognizes handpieces for easy set-up
4. Exceptional control of specifically design tips enabling soft tissue and fine bone dissection in close proximity to delicate anatomic structure
5. Light weight, ergonomic design provides superb tactile feedback and requires no external cooling
6. Suction, Irrigation and Ultrasonic functions can be used alone or simultaneously by using the SYNC mode
7. Large selection of specialty tips available

B. SPECIFICATIONS

1.CONSOLE

Dimensions & Weight

- Width: 33.02 cm
- Height: 23.37 cm
- Depth: 44.20 cm
- Weight: 18.28 kg

Power supply

- Input voltage: 100-240V ~ 50/60 Hz; 6-3A

Fuses

- Type: 2 x 6.3A
- Rating: 250 VAC
- Size: 5 x 20 mm
- Breaking capacity: (F) Fast-Acting (L) Low Breaking Capacity 63 A at 250 VAC, IEC 60127

4 Channel RFID

- Frequency of Operation: 13.56MHz A38
- RF Bandwidth: 13.553 MHz – 13.567 MHz
- Modulation: ASK
- RF Field Strength: 67.92 dBuV/m at 3m

2.Handpiece

Dimensions & Weight

- Housing: 12.4 cm length x 3.0 cm outside diameter
- Cable: 4.6m length
- Connector: 6.9 cm length x 1.8 cm outside diameter
- Weight: 320 g

Frequency

- 25.4 kHz

3. Hand Controller

Housing (L x W x H): 16.8 cm x 5.1 cm x 1.5cm

Cable: 4.6 m length

Connector: 1 cm length x 1.3 cm outside diameter

Weight: 230 g

C. CONSIST OF

1. Ultrasonic Aspirator Console	1ea
2. Sonopet Cart	1ea
3. Advanced Foot pedal	1ea
4. Universal Angled Handpiece	2ea
5. Sterilization Tray	2ea
6. Torque Wrench, 7mm	2ea
7. Hand Controller	1ea

D. Remarks :

The installation and operating should be carried out by the supplier's responsibility at enduser's site.

Three years warranty (Handpiece 1 years) after installation and correct operation.